



受控

INSTRUCTIONS FOR SCALP VEIN SETS

【Scope of Application】

- This product is used together with IV set.

【main structure and properties】

- The product is made up of protective cap, needle, needle stand, tubing and connector;
- Tubing is transparent which can easy to observe;
- The design of needle stand to make it convenient for handle and usage;
- Tubing is flexible enough to against winding & twisting.

【Specification】

- TKSVS-01, TKSVS-02, TKSVS-03, Needle size: 18G to 27G

【Usage】

- Firstly, check valid period and completeness os products before using;;
- Open the package from the peel site and take out the scalp vein set;
- Connecting luer lock connector of needle with infusion set;
- Remove protective cap of needle, Then do as sterile aseptic manipulation and fix it.

【Precautions】

- It may causes hypothermia-esp. Weak patients;
- It may causes abnormal rhythm of ventrickes of heart because of fast flow rate;
- There should be guardian during infusion in case of any accident;
- It should be sterile and free from pyrogen before package opened;
- It has a valid period of five years, Do not use if it is beyond valid date;
- Discard after singe use.

【Warnigs and Indications】

- Do not use if the package is opened or damaged, or protective cap is apated, any part is missed, And foreign matter is found;
- It should be used by special trained person, and the skin should be correctly sterilized.
- The products can not be resued.The repetition use may cause cross infection, even serious illnesses such as hepatitis,HIV and other infectious disease.

【Storage】

- Product should be kept in a condition that relative humidity is less than 80%, it also should be free from corrosive gas, cool, dry with good ventilation and clean.

Unique device identifier	Batch code	Catalogue number	Caution	Sterilized using ethylene oxide
Date of manufacture	Manufacturer	CONTAINS or PRESENCE of PHTHALATES	Do not re-use	Do not use if package is damaged
Medical device	Use-by date	Do not sterilize	Fragile, handle with care	
Keep away from sunlight	Keep dry	Non-pyrogenic	This way up	
Single sterile barrier system	Single sterile barrier system with protective packaging outside		REV.: A/0 DATE: 2023-11-14	