

Instruction Of Extension Tubes

【Scope of Application】

- The product is applied in clinical infusion in the operation,there is need to be connected to infusion containers,a variety of drug infusion at the same time,the extension of the infusion pipeline.

【Main structure and properties】

- Tube is transparent,with normal or corrected vision to find bubbles in line ;
- Tube is flexible enough to resist winding and twisting;
- The design of closure-piercing device to make it convenient for handle and usage;
- All interfaces should be cases,cases should not fall off naturally,and easy to dismantle.
- After the clamps hold the canula,should have no liquid flow,and should not damage the catheter in continuous use.

【Specification】

- TKET-F01, TKET-F02, TKET-F03.

【Usage】

- Firstly, check valid period and completeness of package before using;
- Tear the package at the cut part and take out the extension tubes;
- Extension tubes joint connection on the infusion vessel joint;
- According to the aseptic operation,start infusion or infusion shdrugs.
- Infusion operation must have the professionals nursing to implement and guardian.

【Precautions】

- It may causes hypothermia-esp, weak patients;
- It may causes abnormal thymem of ventricketes of heart because of fast flow rate;
- There should maintain the guardianship during the infusion in case of any accident;
- It should be sterile and free from pyrogen before package opened;.
- It has a valid period of three years, Do not use if it is beyond valid date.
- Discard after singe use.

【warnings and indications】

- Do not use if the package is opened or damaged, or protective cap is aparted, any part is missed, And foreign matter is found;
- It should be used by special trained person.
- The products can not be resued.The repetition use may cause cross infection, even serious illnesses such as hepatitis,HIV and other infectious disease.

【storage】

- Product should be kept in a condition that relative humidity is less than 80%, it also should be free from corrosive gas, cool, dry with good ventilation and clean.

ISO 8536-9-CL-P

Unique device identifier	Batch code	Catalogue number	Caution	Sterilized using ethylene oxide
Date of manufacture	Manufacturer	DOES NOT CONTAINS or PRESENCE of PHTHALATES	Do not re-use	Do not use if package is damaged
Medical device	Use-by date	Do not re-sterilize	Fragile, handle with care	
Keep away from sunlight	Keep dry	Non-pyrogenic	This way up	
Single sterile barrier system		Single sterile barrier system with protective packaging outside		REV.: A/0 DATE: 2023-11-14